

UltraFlux

PCR Plates and Strip Tubes
for Molecular Diagnostic Kits


Designed specifically for use with ABI and Roche Thermal Cyclers

UltraFlux PCR Plates and Strip Tubes have been the choice of kit manufacturers around the world for more than 25 years. Each is designed to meet or exceed the specifications and performance of plates supplied by the original equipment manufacturer.

Advantages of UltraFlux PCR Products

UltraFlux PCR Plates

High-contrast, black alphanumeric indexing.




Raised wells ensure optimal sealing to prevent evaporation.

Ultra clear wells for easy sample visualization.

Barcoding available for select items.



UltraFlux PCR Strip Tubes



UltraFlux triple interwell linkages make the tubes rigid, decreasing the chance of spillage.

White tubes with optically-clear flat caps. White wells enhance fluorescence reflection while optically-clear, flat caps ensure that the fluorescence signal is passed through uninhibited.

Why buy kit consumables directly from the manufacturer?

- Cost savings
- Certified lot testing
- Improved logistics
- Supply chain control



A trusted partner to diagnostic kit manufacturers around the globe.

ISO 8 / Class 100,000 Cleanroom production and packaging ensures purity.



Injection molding expertise and in-house mold craftsmen ensure repeatable performance and reliable supply.



Certified lot testing under a stringent ISO-9001 Quality System guarantees consistent, reliable product.



Source your kit consumables from a trusted leader for improved performance and profitability.



www.SSIbio.com



Made in USA

UltraFlux

PCR Plates and Strip Tubes
for Molecular Diagnostic Kits

Designed specifically for use with ABI and Roche Thermal Cyclers

UltraFlux PCR Plates and Strip Tubes have been the choice of kit manufacturers around the world for more than 25 years. Each is designed to meet or exceed the specifications and performance of plates supplied by the original equipment manufacturer.



ABI Fast® (0.1 mL)



3426-00



3145-00



3248-00

ABI Standard (0.2 mL)



3425-00



3135-00



3247-00



Roche LightCycler®, White



3455-40



3145-40



3248-40

Roche LightCycler®, Clear



3455-00



3145-00



3248-00

UltraFlux PCR

| ABI Fast | |
|--------------|--|
| 3426-00 | 0.1 mL Low-Profile, 96-Well Plate, Clear |
| 3145-00 | 0.1 mL Low-Profile 8-Strip Tubes & Caps, Clear |
| 3248-00 | 0.1 mL Low-Profile 8-Strip with Attached Caps, Clear |
| 3248-40 | 0.1 mL Low-Profile 8-Strip with Attached Caps, White |
| ABI Standard | |
| 3425-00 | 0.2 mL 96-Well Plate, Clear |
| 3135-00 | 0.2 mL 8-Strip Tubes, Clear |
| 3247-00 | 0.2 mL 8-Strip Tubes with Attached Caps, Clear |

All caps are optically clear and flat, making them ideal for qPCR.

Fast is a trademark of Applied Biosystems.
LightCycler is a registered trademark of Roche Diagnostics, GmbH.

| Roche LightCycler, White | |
|--------------------------|--|
| 3455-40 | 0.1 mL Low-Profile, 96-Well Plate, White |
| 3145-40 | 0.1 mL Low-Profile 8-Strip Tubes & Caps, White |
| 3248-40 | 0.1 mL Low-Profile 8-Strip Tubes with Attached Caps, White |
| Roche LightCycler, Clear | |
| 3455-00 | 0.1 mL Low-Profile, 96-Well Plate, Clear |
| 3145-00 | 0.1 mL Low-Profile 8-Strip Tubes & Caps, Clear |
| 3248-00 | 0.1 mL Low-Profile 8-Strip Tubes with Attached Caps, Clear |



Contact me now
for a
FREE
UltraFlux
Evaluation Kit

Steve Laramore
Director of Sales, Americas

Mobile: +1 (614) 947-9591
S.Laramore@SSIbio.com

