



SECTION A-A



| ITEM | PART NUMBER | DESCRIPTION                   | QTY. |
|------|-------------|-------------------------------|------|
| 1    | 6131-XX     | 3.6mL INT. THD. CRYOVIAL TUBE | 1    |
| 2    | M9-0-2285   | SILICONE CAP SEAL             | 1    |
| 3    | 6100-XX     | EXT. THD. CRYOVIAL CAP        | 1    |

|  |        |  |         |           |             |             |
|--|--------|--|---------|-----------|-------------|-------------|
| CAD GENERATED DRAWING.<br>DO NOT MANUALLY UPDATE |        |  <b>1310 Thurman Street</b><br><b>Lodi CA, 95240</b><br><b>Ph: (209) 333-2120</b><br><b>Fx: (209) 333-8623</b><br><b>info@scientificspecialties.com</b> |         |           |             |             |
| APPROVALS  | DATE   |  |         |           |             |             |
| DRAWN<br>S.Gallagher                             | 9/1/11 | <div style="text-align: center;"> <h1>6131-S0</h1> <h2>3.6mL CRYOVIAL ASSEMBLY</h2> </div>   |         |           |             |             |
| CHECKED  |        |  |         |           |             |             |
| RESP ENG   |        |  |         |           |             |             |
| MFG ENG  |        |  |         |           |             |             |
| FINISH   | ---    | QUAL ENG   |         | SIZE<br>A | SHRINK RATE | REV.<br>REL |
| DO NOT SCALE DRAWING                             |        | SCALE  | JOB NO. | SHEET     | 1 OF 1      |             |

PRINTED 9/2/2011