

Macro Volume Tip Compatibility Chart

5ml Tips

| Brand | 4410/4411 | 4420/4421 |
|-----------|-----------|-----------|
| Eppendorf | | X |
| Gilson | X | |
| Rainin | X | |
| Socorex | | X |

10ml Tips

| Brand | 4520/4521 | 4530/4531 |
|-------------|-----------|-----------|
| Eppendorf | | X |
| Gilson | | X |
| Rainin | | X |
| Socorex | X | |
| Accumax | X | |
| Biohit | X | |
| Finnpipette | X | |

5ml Tips

| | | | | |
|----------------|--|--|----------------------------|--|
| 4410-00 | 5ml, Gilson-fit, Bulk, Natural | | B | 250 per Pack, 10 Packs per Case |
| 4420-00 | 5ml, Eppendorf-fit, Bulk, Natural | | B | 250 per Pack, 10 Packs per Case |
| 4411-00 | 5ml, Gilson-fit, Racked, Natural | | R | 50 Tips per Rack, 10 Racks per Case |
| 4421-00 | 5ml, Eppendorf-fit, Racked, Natural | | R | 50 Tips per Rack, 10 Racks per Case |
| 4411-S0 | 5ml, Gilson-fit, Racked, Sterile, Natural | | R S | 50 Tips per Rack, 10 Racks per Case |
| 4421-S0 | 5ml, Eppendorf-fit, Racked, Sterile, Natural | | R S | 50 Tips per Rack, 10 Racks per Case |
| 4411-SF | 5ml, Gilson-fit, Racked, Sterile, Filtered, Natural | | R S F | 50 Tips per Rack, 10 Racks per Case |
| 4421-SF | 5ml, Eppendorf-fit, Racked, Sterile, Filtered, Natural | | R S F | 50 Tips per Rack, 10 Racks per Case |

10ml Tips

| | | | | |
|----------------|--|--|----------------------------|--|
| 4520-01 | 10ml, Bulk, Natural | | B | 200 per Pack, 5 Packs per Case |
| 4530-00 | 10ml, Gilson-fit, Bulk, Natural | | B | 200 per Pack, 5 Packs per Case |
| 4521-01 | 10ml, Racked, Natural | | R | 25 Tips per Rack, 10 Racks per Case |
| 4531-00 | 10ml, Gilson-fit, Racked, Natural | | R | 25 Tips per Rack, 10 Racks per Case |
| 4521-S1 | 10ml, Racked, Sterile, Natural | | R S | 25 Tips per Rack, 10 Racks per Case |
| 4531-S0 | 10ml, Gilson-fit, Racked, Sterile, Natural | | R S | 25 Tips per Rack, 10 Racks per Case |
| 4521-SF | 10ml, Racked, Sterile, Filtered, Natural | | R S F | 25 Tips per Rack, 10 Racks per Case |
| 4531-SF | 10ml, Gilson-fit, Racked, Sterile, Filtered, Natural | | R S F | 25 Tips per Rack, 10 Racks per Case |