

Pipettor Compatibility Chart

NON-STERILE AND STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117N00S (4117NS0S)	4137N00S (4137NS0S)	4237N00S (4237NS0S)	4287N00S (4287NS0S)	4927N00S (4927NS0S)	4337N00S (4337NS0S)	4347N00S (4347NS0S)	
Brand	Transferpette S	1 µL	0.1 - 1 µL	SINGLE	4117NS0S	4137NS0S						
		2.5 µL	0.1 - 2.5 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S					
		100 µL	10 - 100 µL	SINGLE			4237NS0S					
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
		200 µL	20 - 200 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
		300 µL	15 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
Capp	Aero	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S			
		50 µL	5 - 10 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	25 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		300 µL	50 - 300 µL	SINGLE					4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
			10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S					
			50 µL	5 - 50 µL	MULTI			4237NS0S	4287NS0S	4927NS0S		
	Bravo	2 µL	0.1-2 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S			
		50 µL	5 - 50 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
1000 µL		100 - 1000 µL	SINGLE							4337NS0S	4347NS0S	
			10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S					
		50 µL	5 - 10 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
Eppendorf	Reference	2.5 µL	0.1 - 2.5 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S					
		100 µL	10 - 100 µL	SINGLE			4237NS0S					
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
	Reference 2	2.5 µL	0.1 - 2.5 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
				20 µL	2 - 20 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S	

Pipettor Compatibility Chart

NON-STERILE AND STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117N00S (4117NS0S)	4137N00S (4137NS0S)	4237N00S (4237NS0S)	4287N00S (4287NS0S)	4927N00S (4927NS0S)	4337N00S (4337NS0S)	4347N00S (4347NS0S)	
Eppendorf	Reference 2	100 µL	10 - 100 µL	SINGLE			4237NS0S		4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		300 µL	30 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
	Research	2.5 µL	0.1 - 2.5 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE								
		100 µL	10 - 100 µL	SINGLE			4237NS0S					
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		300 µL	30 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
		10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S						
		100 µL	10 - 100 µL	MULTI			4237NS0S					
	300 µL	30 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S				
	Research plus	2.5 µL	0.1 - 2.5 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL (yellow)	SINGLE			4237NS0S					
		20 µL	2 - 20 µL (light gray)	SINGLE	4117NS0S	4137NS0S						
		100 µL	10 - 100 µL	SINGLE			4237NS0S					
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		300 µL	30 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
		10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S						
	300 µL	20 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S				
	Xplorer / Xplorer Plus	10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	1 - 20 µL	SINGLE			4237NS0S	4287NS0S				
		100 µL	5 - 100 µL	SINGLE			4237NS0S	4287NS0S				
		200 µL	10 - 200 µL	SINGLE			4237NS0S	4287NS0S				
		300 µL	15 - 300 µL	SINGLE					4927NS0S			
		1000 µL	50 - 1000 µL	SINGLE							4337NS0S	4347NS0S
		10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S						
		100 µL	5 - 100 µL	MULTI			4237NS0S	4287NS0S				
		300 µL	15 - 300µL	MULTI					4927NS0S			
1200 µL	50 - 1200µL	MULTI								4347NS0S		
Finnpipette	4500	10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S			

Pipettor Compatibility Chart

NON-STERILE AND STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117N00S (4117NS0S)	4137N00S (4137NS0S)	4237N00S (4237NS0S)	4287N00S (4287NS0S)	4927N00S (4927NS0S)	4337N00S (4337NS0S)	4347N00S (4347NS0S)	
Finnpipette	4500	50 µL	5 - 50 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	200 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
		50 µL	5 - 50 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
	F1	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		5 µL	0.5 - 5 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	1 - 10 µL (pink)	SINGLE	4117NS0S	4137NS0S						
		10 µL	1 - 10 µL (yellow)	SINGLE			4237NS0S		4927NS0S			
		20 µL	2 - 20 µL (micro)	SINGLE	4117NS0S	4137NS0S						
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
	F2	300 µL	30 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		5 µL	0.5 - 5 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	1 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S			
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
	Focus	100 µL	10 - 100 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
		300 µL	30 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
		50 µL	5 - 50 µL	SINGLE			4237NS0S		4927NS0S			
	Novus	300 µL	30 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
		10 µL	1 - 10 µL	MULTI	4117NS0S	4137NS0S						
		50 µL	5 - 50 µL (yellow)	MULTI			4237NS0S		4927NS0S			
		300 µL	30 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
		1200 µL	100 - 1200 µL	MULTI						4337NS0S	4347NS0S	
Gilson		Pipetman Classic	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S					
			10 µL	1 - 10 µL	SINGLE	4117NS0S	4137NS0S					
			20 µL	2 - 20 µL	SINGLE			4237NS0S				
			1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S

Pipettor Compatibility Chart

NON-STERILE AND STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117N00S (4117NS0S)	4137N00S (4137NS0S)	4237N00S (4237NS0S)	4287N00S (4287NS0S)	4927N00S (4927NS0S)	4337N00S (4337NS0S)	4347N00S (4347NS0S)	
Gilson	Pipetman Classic	100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
		20 µL	2 - 20 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
		300 µL	30 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
	Pipetman L/M/G	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	1 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S			
		100 µL	10 - 100 µL	SINGLE			4237NS0S		4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
		10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	MULTI			4237NS0S		4927NS0S			
		200 µL	20 - 200 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
	Pipetman Neo	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	1 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S					
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
1000 µL		100 - 1000 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S				
Pipetman Ultra	10 µL	1 - 10 µL	SINGLE		4137NS0S							
	20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S				
	100 µL	20 - 100 µL	SINGLE			4237NS0S		4927NS0S				
	300 µL	20 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S				
	1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S		
Hamilton	Softgrip	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	1 - 10 µL	SINGLE		4137NS0S						
		25 µL	2.5 - 25 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		300 µL	30 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
		300 µL	30 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
Hirschmann	labopette	3 µL	0.1 - 3 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						

Pipettor Compatibility Chart

NON-STERILE AND STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL
					4117N00S (4117NS0S)	4137N00S (4137NS0S)	4237N00S (4237NS0S)	4287N00S (4287NS0S)	4927N00S (4927NS0S)	4337N00S (4337NS0S)	4347N00S (4347NS0S)
Hirschmann	labopette	20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S		
		100 µL	10 - 100 µL	SINGLE			4237NS0S		4927NS0S		
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S		
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S
HTL Pipette Discovery	Pro / Comfort / LABMATE pro / Optipette	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S					
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S					
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S		
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S		
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S		
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S
IKA PETTE	vario	2 µL	0.1 - 2.0 µL	SINGLE	4117NS0S	4137NS0S					
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S					
		20 µL	2 - 20 µL	SINGLE			4237NS0S				
		100 µL	10 - 100 µL	SINGLE			4237NS0S				
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S		
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S
		10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S					
		200 µL	20 - 200 µL	MULTI			4237NS0S	4287NS0S	4927NS0S		
Nichiryo	Oxford Benchmate / Benchmate II	2 µL	0.1 - 2 µL	SINGLE	4117NS0S	4137NS0S					
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S					
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S		
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S		
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S		
		1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S
	EX II MULTI	10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S					
		100 µL	5 - 100 µL	MULTI			4237NS0S		4927NS0S		
		200 µL	20 - 200 µL	MULTI			4237NS0S		4927NS0S		
		300 µL	30 - 300 µL	MULTI				4287NS0S	4927NS0S		
	Nichipet EX / EX II / EX PLUS II / LE / Premium / Premium LT	10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S					
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S		
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S		
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S		
1000 µL		100 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
Rainin	Pipet-Lite XLS (Universal)	2 µL	0.1 - 2 µL	SINGLE	4117NS0S	4137NS0S					
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S					

Pipettor Compatibility Chart

NON-STERILE AND STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117N00S (4117NS0S)	4137N00S (4137NS0S)	4237N00S (4237NS0S)	4287N00S (4287NS0S)	4927N00S (4927NS0S)	4337N00S (4337NS0S)	4347N00S (4347NS0S)	
Rainin	Pipet-Lite XLS (Universal)	20 µL	2 - 20µL	SINGLE			4237NS0S		4927NS0S			
		100 µL	10 - 100 µL	SINGLE			4237NS0S		4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		300 µL	20 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
Sartorius	(Biohit) mLINE	3 µL	0.1 - 3 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5-10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S			
		100 µL	5 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 -1000 µL	SINGLE							4337NS0S	4347NS0S
		10 µL	0.5 -10 µL	MULTI	4117NS0S	4137NS0S						
		100 µL	5 - 100 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
	300 µL	30 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S				
	(Biohit) Proline Plus / Proline	3 µL	0.1 - 3 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S			
		50 µL	5 - 50 µL	SINGLE			4237NS0S					
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S				
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
		10 µL	0.5 -10 µL	MULTI	4117NS0S	4137NS0S						
		100 µL	10 - 100 µL	MULTI			4237NS0S	4287NS0S				
		300 µL	30 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
	Picus	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.2 - 10 µL	SINGLE	4117NS0S	4137NS0S						
120 µL		5 - 120 µL	SINGLE				4287NS0S	4927NS0S				
300 µL		10 - 300 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S				
1000 µL		50 - 1000 µL	SINGLE							4337NS0S	4347NS0S	
10 µL		0.2 - 10 µL	MULTI	4117NS0S	4137NS0S							
120 µL		5 - 120 µL	MULTI				4287NS0S	4927NS0S				
300 µL		10 -300 µL	MULTI	4117NS0S	4137NS0S							
300 µL		10 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S				
1200 µL		50 - 1200 µL	MULTI							4337NS0S	4347NS0S	
Tacta	3 µL	0.1 - 3 µL	SINGLE	4117NS0S	4137NS0S							
	10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S							

Pipettor Compatibility Chart

NON-STERILE AND STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117N00S (4117NS0S)	4137N00S (4137NS0S)	4237N00S (4237NS0S)	4287N00S (4287NS0S)	4927N00S (4927NS0S)	4337N00S (4337NS0S)	4347N00S (4347NS0S)	
Sartorius	Tacta	20 µL	2 - 20 µL	SINGLE			4237NS0S	4287NS0S				
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S				
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S				
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
		3 µL	0.1 - 3 µL	MULTI	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S						
		100 µL	5 - 100 µL	MULTI			4237NS0S	4287NS0S				
		300 µL	30 - 300 µL	MULTI			4237NS0S	4287NS0S				
Socorex	Acura	2 µL	0.1 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	1 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S					
		50 µL	5 - 50 µL	SINGLE			4237NS0S					
		100 µL	10 - 100 µL	SINGLE			4237NS0S		4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
		10 µL	0.5 - 10 µL	MULTI	4117NS0S	4137NS0S						
		10 µL	1 - 10 µL	MULTI	4117NS0S	4137NS0S						
		50 µL	5 - 50 µL	MULTI			4237NS0S					
		100 µL	10 - 100 µL	MULTI			4237NS0S		4927NS0S			
		200 µL	20 - 200 µL	MULTI			4237NS0S	4287NS0S				
	350 µL	40 - 350 µL	MULTI			4237NS0S	4287NS0S	4927NS0S				
		Calibra 822	2 µL	0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S					
	10 µL		1 - 10 µL	SINGLE	4117NS0S	4137NS0S						
	10 µL		1 - 10 µL (200µL tips)	SINGLE			4237NS0S					
100 µL	10 - 100 µL		SINGLE			4237NS0S	4287NS0S					
200 µL	20 - 200 µL		SINGLE			4237NS0S	4287NS0S					
	Calibra 852	1000 µL	100 - 1000 µL	SINGLE						4337NS0S	4347NS0S	
10 µL		1 - 10 µL	MULTI	4117NS0S	4137NS0S							
100 µL		10 - 100 µL	MULTI			4237NS0S	4287NS0S					
	Ergonomic	200 µL	20 - 200 µL	MULTI			4237NS0S	4287NS0S				
2 µL		0.2 - 2 µL	SINGLE	4117NS0S	4137NS0S							
VWR	Ergonomic	10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S					
		50 µL	5 - 50 µL	SINGLE			4237NS0S	4287NS0S				
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S				
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S				

Pipettor Compatibility Chart

NON-STERILE AND STERILE PIPETTE TIPS



Brand	Model	Max Vol	Volume Range	Channel	10 µL	10 µL XL	200 µL	200 µL XL	300 µL	1000 µL	1250 µL	
					4117N00S (4117NS0S)	4137N00S (4137NS0S)	4237N00S (4237NS0S)	4287N00S (4287NS0S)	4927N00S (4927NS0S)	4337N00S (4337NS0S)	4347N00S (4347NS0S)	
VWR	Ergonomic	200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S				
		250 µL	50 - 250 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S
		10 µL	1 - 10 µL	MULTI	4117NS0S	4137NS0S						
		50 µL	5 - 50 µL	MULTI			4237NS0S	4287NS0S				
		200 µL	20 - 200 µL	MULTI			4237NS0S	4287NS0S				
		300 µL	50 - 300 µL	MULTI			4237NS0S	4287NS0S	4927NS0S			
	Ultra	2 µL	0.1 - 2 µL	SINGLE	4117NS0S	4137NS0S						
		10 µL	0.5 - 10 µL	SINGLE	4117NS0S	4137NS0S						
		20 µL	2 - 20 µL	SINGLE			4237NS0S		4927NS0S			
		100 µL	10 - 100 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		200 µL	20 - 200 µL	SINGLE			4237NS0S	4287NS0S	4927NS0S			
		1000 µL	100 - 1000 µL	SINGLE							4337NS0S	4347NS0S