



SECTION A-A



ITEM	PART NUMBER	DESCRIPTION	QTY.
1	6222-XX	1.8mL EXT. THD. CRYOVIAL TUBE	1
2	6200-XX	INT. THD. CRYOVIAL CAP	1

CAD GENERATED DRAWING. DO NOT MANUALLY UPDATE		 1310 Thurman Street Lodi CA, 95240 Ph: (209) 333-2120 Fx: (209) 333-8623 info@scientificspecialties.com		
APPROVALS	DATE			
DRAWN S.Gallagher	9/1/11	<h1>6222-S0</h1> <h2>1.8mL CRYOVIAL ASSEMBLY</h2>		
CHECKED				
SSI MATERIAL: ---	RESP ENG			
MANUFACTURERS MATERIAL: ---	MFG ENG			
FINISH ---	QUAL ENG	SIZE A	SHRINK RATE	REV. REL
DO NOT SCALE DRAWING		SCALE	JOB NO.	SHEET 1 OF 1