





NOTE(S): UNLESS OTHERWISE SPECIFIED.

1. MATERIAL: 2.0mil CLEAR POLYESTER FILM WITH 1.0-1.3mil ACRYLIC ADHESIVE ON A 2.0mil CLEAR POLYESTER RELEASE LINER.
2. BACKSCORES THROUGH 2.0mil CLEAR POLYESTER RELEASE LINER ONLY.
3. PERFORATIONS THROUGH 2.0mil POLYESTER FILM ONLY.



UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETERS	CAD GENERATED DRAWING, DO NOT MANUALLY UPDATE		 1310 Thurman Street Lodi CA, 95240 Ph: (209) 333-2120 Fx: (209) 333-8623 info@scientificspecialties.com
	APPROVALS	DATE	
SSI MATERIAL:	DRAWN: P. DO	2/21/2018	3612-00-XL THERMALSEAL RT2
MANUFACTURER'S MATERIAL:	CHECKED:		
FINISH:	DESIGN ENG:		SIZE: SHRINK RATE: A
DO NOT SCALE DRAWING	MFG ENG:		
	QUAL ENG:		SCALE: 1:2 JOB NO.: SHEET: 1 OF 1
			REV.: REL

PRINTED 2/22/2018