



**SECTION A-A**  
**SCALE: 2.2 : 1**

UNLESS OTHERWISE SPECIFIED  
DIMENSIONS ARE IN MILLIMETERS  
TOLERANCES ARE:  
DECIMALS ANGLES  
.XX ± .25 ± 0.1  
.XXX ± .025  
.XXXX ± .0127

CAD GENERATED DRAWING.  
DO NOT MANUALLY UPDATE

APPROVALS	DATE
DRAWN: P. DO	3/10/2015

SSI MATERIAL:  
----

CHECKED:

MANUFACTURER'S MATERIAL:  
----

RESP ENG:

FINISH:  
----

MFG ENG:

DO NOT SCALE DRAWING

QUAL ENG:



1310 Thurman Street  
Lodi CA, 95240  
Ph: (209) 333-2120  
F: (209) 333-8623  
info@scientificspecialties.com

**1310-00 (METRIC)**  
2.0 mL "GRADUATED" TUBE W/ CAP

SIZE: A	SHRINK RATE: ----	REV.: 3
SCALE: 1:1	JOB NO.:	SHEET: 1 OF 1

PRINTED 3/11/2015